

iFusion

- MATLAB prototype.
- Two layers:
 - INTERFACE: Dedicated to Human Computer Interaction routines (11 functions)
 - CORE: Dedicated to the REAL part of the iFusion project (12 functions)
 - IMAge Processing
 - GEOmetry
 - MISCellaneous

Human Computer Interface (HCI)

Ifusion_main.m

Ifusion_Global.m

Ifusion_interface_1.m

Ifusion_interface_2.m

Ifusion_interface_3.m

Ifusion_CF_box.m

Ifusion_Continue_box.m

Ifusion_Warning_box.m

Ifusion_Control_Points.
m

Ifusion_Trace_Curves.m

UserInteractionCurve.m

catetherPathFast.m

detectCatheterTipMultipleFrames.m

imageEnhancement.m

solveEq18.m

IMage Processing

Project_point_line.m

Epipolar.m

StretchLine.m

Intersection_plane_line.m

angles2refsys.m

Intersection_line_line.m

ComputeAngles.m

GEOMETRY

MISCellaneous

extractEndDiastole.m

CORE functions

iFusion: INTERFACE

iFusion: INTERFACE

- Similar to the one of the first C++ prototype.
- Actual visualization (buttons, menus, etc.) quite differs from the previous version.
- NO functions or methods regarding the iFusion project.

iFusion: INTERFACE_1

Loading Angiographies

Ifusion_interfaz_1

Load Case

or SELECT THE ANGIOGRAPHY SEQUENCES:

BPB (Left)	BPB (Right)	APB (Right)	APB (Left)	Coronary
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Load	Load	Load	Load	Load

OK

iFusion: INTERFACE_1

Loading Angiographies

Ifusion_interfaz_1

Load Case

or SELECT THE ANGIOGRAPHY SEQUENCES:

BPB (Left)	BPB (Right)	APB (Right)	APB (Left)	Coronary
Load	Load	Load	Load	Load

Previews of loaded Angiographies

OK

iFusion: INTERFACE_1

Loading Angiographies

Ifusion_interfaz_1

Load Case

or SELECT THE ANGIOGRAPHY SEQUENCES:


An entire clinical case
can be Loaded

BPB (Left)	BPB (Right)	APB (Right)	APB (Left)	Coronary
Load	Load	Load	Load	Load

OK


iFusion: INTERFACE_2

Curve Model Generation




iFusion: INTERFACE_2

Curve Model Generation




iFusion: INTERFACE_2

Curve Model Generation




iFusion: INTERFACE_2

Curve Model Generation




iFusion: INTERFACE_2

Curve Model Generation




iFusion: INTERFACE_2

Curve Model Generation




iFusion: INTERFACE_2

Curve Model Generation




iFusion: INTERFACE_2

Curve Model Generation




iFusion: INTERFACE_2

Curve Model Generation




iFusion: INTERFACE_2

Curve Model Generation




iFusion: INTERFACE_2


Curve Model Generation




iFusion: 3D vessel reconstruction




iFusion: complete 3D geometry of the global and local reference systems



iFusion: complete 3D geometry of the global and local reference systems




IVUS to Angio correspondence




iFusion: INTERFACE_3

IVUS to Angio correspondence




iFusion: 3D visualization of vessel and selected positions



iFusion: INTERFACE

- *Ifusion_main*: Define Global variables and call the three interfaces in the correct order.
- *Ifusion_Global*: Contains global variables used by all the INTERFACE functions.

Human Computer Interface (HCI)



iFusion: INTERFACE

- *Ifusion_interface_1*: Load the X-ray DICOM files.
- *Ifusion_interface_2*: Interact with the user to find starting and ending points of the desired structure to be modeled; generate 3D model.
- *Ifusion_interface_3*: Load the pullback and visualize results.

Human Computer Interface (HCI)

Ifusion_main.m

Ifusion_Global.m

Ifusion_interface_1.m

Ifusion_interface_2.m

Ifusion_interface_3.m

Ifusion_CF_box.m

Ifusion_Continue_box.m

Ifusion_Warning_box.m

Ifusion_Control_Points.
m


Ifusion_Trace_Curves.m

UserInteractionCurve.m

iFusion: INTERFACE

- *Ifusion_CF_box*: Interaction with the user to input the CF distance if not provided by the infos of the DICOM files.
- *Ifusion_Continue_box*: Yes/No box to assure the user desires to continue with the processing.
- *Ifusion_Warning_box*: Warning message if DICOM files are missing.

Human Computer Interface (HCI)



iFusion: INTERFACE

- *Ifusion_Control_Points*: Given two paths (left and right) appropriately assigns Control Points.
- *Ifusion_Trace_Curves*: Draw model curves on the display.
- *UserInteractionCurve*: Update the curves when a control point is moved.

Human Computer Interface (HCI)

Ifusion_main.m

Ifusion_Global.m

Ifusion_interface_1.m

Ifusion_interface_2.m

Ifusion_interface_3.m

Ifusion_CF_box.m

Ifusion_Continue_box.m

Ifusion_Warning_box.m

Ifusion_Control_Points.
m

Ifusion_Trace_Curves.m

UserInteractionCurve.m

iFusion: CORE

iFusion: CORE

- IMP (Image Processing): 4 functions.
- GEO (Geometry): 7 functions.
- MISC (Miscellaneous): 1 function.

iFusion: CORE, IMP

- *catheterPathFast*: Implements the Global Minimum for Active Contour Models (see Cohen et. Al.).
- *detectCatheterTipMultipleFrames*: Detects rounded objects in a sequence of frames.
- *imageEnhancement*: Enhances the image by local histogram equalization.



- *solveEq18*: Specific routine used by the method *catheterPathFast*. Actually rapidly solves Eq. 18 in Cohen et Al.

iFusion: CORE, GEO

- *Project_point_line*: Projects a point to the closest point of a line, in 2D.
- *Epipolar*: Builds the epipolar line of a reference system, with respect to another.
- *StretchLine*: Takes an epipolar segment and crops it to fit with the image size.



iFusion: CORE, GEO

- *Intersection_plane_line*: Finds the 3D point that results from the intersection of a given plane and a line.
- *angles2refsys*: Computes the local reference system starting from CArm angles.
- *Intersection_line_line*: Finds the 3D point that is the closest of both given lines, or if it exists, the intersection between lines.
- *ComputeAngles*: Computes the angle of a 2D vector, i.e. the angle of polar coordinates of the vector.




iFusion: CORE, MISC

- *extractEndDiastole*: If the ECG signal is available, this function is used to perform an ECG gating at end diastole. It can be easily modified to perform the gating in different ECG positions.


iFusion: Calls architecture

Schema of calls: ALL



Schema of calls: Intra-blocks

Human Computer Interface (HCI)



detectCatheterTipMultipleFrames.m

catetherPathFast.m

imageEnhancment.m

solveEq18.m

IMage Processing

CORE functions

StretchLine.m

angles2refsys.m

GEOMETRY

Epipolar.m

Intersection_line_line.m

Intersection_plane_line.m


Project_point_line.m

ComputeAngles.m

MISCellaneous

extractEndDiastole.m

Schema of calls: Inter-blocks in RED




iFusion: Flow diagrams and Graphical explanations

(only of some functions)


Flow Diagram

(HCI: UserInteractionCurve.m)




Flow Diagram

(HCl: Ifusion_Control_Points.m)




Graphical explanations (HCl: StretchLine.m)



Graphical explanations


(HCl: Intersection_line_line.m)





Graphical explanations

(Find optimal point pair)

Left Angiography



Right Angiography



3D anti-projection error
(see 'e' in previous slide)

Best pair:
- Left: 8
- Right: 5

iFusion: In Code Documentation

(very detailed for all CORE functions)

(sufficient for all HCI functions)

Function Header Comments

- SYNOPSIS
- INPUT variables
- OUTPUT variables
- REFERENCES
- COMMENTS

Example

```
>> help imageEnhancement
```

```
SYNOPSIS: Enhances the input image by means of local histogram  
          qualization on a 27x27 neighborhood.
```

```
INPUT:    im_in:  an image
```

```
OUTPUT:    im_out: enhanced image
```

```
REF:       - Fusion.pdf - Page 6, Sect 2.1  
            - Fusion.pdf - Page 7, Figure 2.3  
            - Equalization as in "MULTIMODAL REGISTRATION OF INTRAVASCULAR  
              ULTRASOUND IMAGES AND ANGIOGRAPHY", page 2, formula (1).
```

```
COMMENTS: This algorithm can be performed by means of the following  
          (commented) for loop. However, the matlab function nlfilter  
          allows to do it a bit faster.
```

Functions add-ons

- Flag 'visualize' (normally set to FALSE)
- If visualize = TRUE
 - open one or more figures
 - visualize the algorithm computation graphically.

iFusion: Technical requirements

Technical requirements and limitations

- The CODE has been developed on:
 - MATLAB 7.7.0.471 (R2008b)
 - Microsoft XP Version 5.1 (Build 2600: Service Pack 3)
 - Java VM Version: 1.6.0_04 with Sun Microsystems Inc. Java HotSpot™ Client VM.
- If the PullBack is too big (causing an Out-of-Memory error), the third interface works in 'simulation' mode, thus not visualizing the pullback.

iFusion: Credits

Original software developer: *David Rotger, PhD*

Leading Technology transfer: *Juan Diego Gomez, Ph.D. Student. and Carlo Gatta, Ph.D.*

MATLAB Coding:

Juan Diego Gomez, Ph.D. Student.

Head of the group: *Petia Radeva Ivanova, Ph.D.*